

*

SPECIFICATION OF DEXTRAN G-25

Items	Descriptions
Appearance	white, slightly off white, dry powder
Observation	no visual color particles
Odor	no obvious smell
Marix	Dextran, cross-linked, formation=beads
Particle size	50~150 μm , (dry bead size)
Swelling rate	x2.5 (2.5 times of powder weight)
Major application	gel filtration
Fractionation range	1000-5000

Shanghai Ren Kai Biotech Co. Ltd.

R-2663, Tenfei Building

#333 Tian Yue Qiao

Shanghai, China

Fax: 8621 6470 2324

Email: novelab@yahoo.com

www.novelab.com

www.ProteinPurify.com